
Validation of cough monitoring by Albus Home, a contactless bedside device for nocturnal monitoring

W. Do¹, R. Russell², C. Wheeler³, H. Javed¹, C. Dogan¹, G. Cunningham¹, V. Khanna¹, M. De Vos⁴, I. Satia⁵, M. Bafadhel⁶, I. Pavord² ¹Albus Health, UK, ²Respiratory Medicine Unit, University of Oxford, UK, ³Royal Brompton Hospital, UK, ⁴Department of Electrical Engineering, KU Leuven, Belgium, ⁵Department of Nedicine, McMaster University, Canada⁶School of Immunology and Microbial Sciences, King's College London, UK.

- Introduction and Aims Methods Adults with respiratory conditions underwent overnight monitoring using Albus Home in their usual bedroom environments. • Objective cough frequency is a key clinical endpoint but existing Participants set-up the plug-and-play device themselves. wearable monitors are limited to 24-hour recordings, which may For reference counts, each audio recording was counted by two annotators, and cough defined as explosive phases be impacted by day-to-day variability. audio-visually labelled by both. • Albus Home uses **contactless** motion, acoustic and • In parallel, recordings were processed by a proprietary Albus system, comprising a deep-learning algorithm with a human environmental sensors to monitor **multiple metrics**, including screening step for verifying coughs or excluding occasional events that mimic cough. respiratory rate and cough, at night-time for as long as required without burdening patients. Performance of the Albus system in detecting individual cough events and reporting hourly cough counts was compared against reference counts. • The aim of this study was to evaluate the performance of Albus Home compared to human manual counts for detecting Accuracy results were reported as sensitivity, specificity, positive predictive value (PPV), negative predictive value (NPV). spontaneous coughs in subjects with respiratory conditions. • In addition, agreement in hourly cough frequencies were analysed using Bland-Altman analysis and Intra-class Correlation Coefficient between Albus system and reference counts.

- observers and reference counts (ICC 0.98, 0.99 respectively).



European Respiratory Society (ERS) International Congress 2022 Disclosure: This research was funded by Albus Health, a spinout from Oxford University, for whom the presenting author is medical director.

Fig. 1. Example use-case and placement of Albus Home RD

Results


30 nights from 10 subjects (COPD=5, asthma=3, sarcoidosis=1, CF=1) comprised 375 hours(h) of recording. Albus counts cough frequencies, with median sensitivity, specificity, PPV and NPV of 94.8, 100.0, 99.1 and 100.0% respectively. • Agreement between Albus and reference was strong (Intra-class Correlation Coefficient 0.99; 95%CI 0.99-0.99; p<0.001) and equivalent to agreement between

s were accurate across l	hours with	high and I	OW
--------------------------	------------	------------	----

ALBUS

HEALTH

Conclusions



• Albus Home provides a **contactless** and **accurate** system for nocturnal cough monitoring, with performance evaluated in an **unseen validation set** of subjects with respiratory conditions, monitored in real-world conditions.

• With the additional potential to concurrently capture a variety nocturnal metrics, the multi-sensory Albus system enables nightly collection of high quality data with applications across clinical care and research.











